

Flash conditions

Device: puriFlash® 4250-iELSD (or now puriFlash® 5.250 pack iELSD)

Solvents: A: Water

B: Acetonitrile

Column: PF-15C18HQ-F0004

Flow rate: 5mL/min

Injection mode: Liquid injection


Crude sample: 50mg

Detection: UV Scan 200-600nm (blue),
ELSD (35°C, Automatic gain) (black)

Pressure: 3bar

Elution conditions:

t (min)	A (%)	B (%)
00:00	95	5
02:00	95	5
05:00	70	30



To achieve this purification:

You will need

- puriFlash® 5.250
[Discover it](#) [Add to card](#)
- Integrated ELSD
[Discover it](#) [Add to card](#)
- puriFlash® column PF-15C18HQ-F0004
[Discover it](#) [Add to card](#)

We highly recommend

- 1mL Stainless Steel Loop with RFID tag
 AYHETO [Add to card](#)
- Ballasting kit for 1/8" tubing - 5 units
 DZ7360 [Add to card](#)
- Manometer ELSD
 FJ6720 [Add to card](#)

