


1. TLC method development



Mobile phase:
86 % Dichloromethane / MeOH 14 %

Compound of interest:
Compound 2

| Compound | Rf | CV |
|----------|------|------|
| 1 | 0.9 | 1.11 |
| 2 | 0.63 | 1.59 |
| 3 | 0.01 | 100 |

$$\Delta CV_{2-1} = 0.48$$

2. Genius proposal



Among the columns proposed by Genius, we have selected a PF-15SIHP-F0012, which was available in stock :

TLC to Flash & Prep (Normal Phase)

Crude sample mg
 Solvent 1 % Additive
 Solvent 2 % Additive



Select your compound(s) of interest

Rf 0.9

Rf 0.63

Rf 0.01

SELECT A COLUMN

| Column | Stock |
|-----------------|--------------------------------------|
| PF-30SIHP-F0025 | ● |
| PF-15SIHP-F0012 | ● |
| PF-30SIHP-F0040 | ● |
| PF-30SIHP-F0080 | ● |
| PF-30SIHP-F0120 | ● |
| PF-30SIHP-F0220 | ● |

Top 3 columns SIHC : for MW < ...

INJECTION MODE

The crude is fully soluble in the below conditions:

Dichloromethane 98.5%
Methanol 1.5%

Yes No

Liquid mL (max 1.66 mL)

Dry load

3. Flash conditions

Device: puriFlash® 5.250P

Solvents: A: Dichloromethane
B: Methanol

Column: PF-15SIHP-F0012

Flow rate: 15mL/min

Injection mode: Liquid injection

Injection volume: 1.5mL

Concentration: 66.67 mg/mL

Crude sample: 100mg


Detection: UV 220nm (purple)

UV Scan 200-400nm (green)

Pressure: 5bar

Elution conditions:

| CV | A (%) | B (%) |
|-------|-------|-------|
| 0.00 | 98.5 | 1.5 |
| 2.50 | 98.2 | 1.8 |
| 6.00 | 96.4 | 3.6 |
| 11.50 | 83.7 | 16.3 |
| 14.50 | 67.6 | 32.4 |
| 17.34 | 67.6 | 32.4 |
| 17.56 | 59 | 41 |
| 18.00 | 0 | 100 |
| 22.00 | 0 | 100 |



To achieve this purification:

You will need

- puriFlash® 5.250P
[Discover it](#) [Add to card](#)
- puriFlash® column PF-15SIHP-F0012
[Discover it](#) [Add to card](#)

We highly recommend

- Safety solvent caps kit - 4 units
BODANO [Add to card](#)
- Safety waste cap with container 5L + Filter
B1SUJ0 [Add to card](#)
- Smartphone for TLC to Flash & Prep application
PHONE0 [Add to card](#)

Download our App

"TLC to Flash & Prep Chromatography"
to make your TLC developments easier
and faster.

